



Product specification

Tube segment opener, for tubes with inner Ø 8.2 mm, HD-PE, white



Product information

Product number	95.1000.800
----------------	-------------

Product descriptions

Name	Tube segment opener, for tubes with inner Ø 8.2 mm, HD-PE, white
Product number	95.1000.800
Product description	Tube segment opener, for tubes with inner Ø 8.2 mm, material: HD-PE, white, 500 piece(s)/case

Material & colours

Color of product	white
------------------	-------

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	216000

Height of inner box	129 mm
Width of inner box	129 mm
Depth of inner box	191 mm
Height of case	150 mm
Width of case	268 mm
Depth of case	398 mm
Volume of case	0.016 cbm
Weight of case	3.6 kg
EAN inner box	4038917124645
EAN case	4038917053273

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT