



## Product specification

### Screw cap micro tube, 0.5 ml, sterile



## Product information

Product number	72.730.027
----------------	------------

## Product descriptions

Name	Screw cap micro tube, 0.5 ml, sterile
Product number	72.730.027
Product description	Screw cap micro tube, working volume: 0.5 ml, skirted conical base, with knurling, transparent, cap: blue, cap assembled, sterile, 100 piece(s)/bag

## Material & colours

Color of product	transparent
O-Ring material	Ethylene propylene diene rubber (EPDM)
Closure material	Polypropylene (PP)
Color of cap	blue

## Size

Diameter	10.8 mm
Length of product	44 mm

## Product characteristics

Base shape	skirted conical base
Graduation	no

Knurls	yes
Cap	standard screw cap
Closure type	Screw cap
Centrifugation max (RCF)	20000 x g

## Purity & certification

Satisfies the requirement	IATA
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	96000
Height of case	140 mm
Width of case	278 mm
Depth of case	348 mm
Volume of case	0.0135 cbm
Weight of case	1.79 kg
EAN inner box	4038917119405
EAN case	4038917043250

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT