



Product specification

Tube segment opener, mounted on tube 75 x 12 mm, HD-PE, white, tube: PS



Product information

Product number	95.1000.476
----------------	-------------

Product descriptions

Name	Tube segment opener, mounted on tube 75 x 12 mm, HD-PE, white, tube: PS
Product number	95.1000.476
Product description	Tube segment opener, mounted on tube 75 x 12 mm, material: HD-PE, white, tube: PS, 50 piece(s)/StackPack

Material & colours

Color of product	white
Tube material	Polystyrene (PS)

Size

Diameter	12 mm
Length of product	75 mm

Product characteristics

Design	mounted on tube 75 x 12 mm
Base shape	Round base

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	400
Piece(s) / pallet	57600
Height of inner box	64 mm
Width of inner box	96 mm
Depth of inner box	146 mm
Height of case	145 mm
Width of case	213 mm
Depth of case	323 mm
Volume of case	0.01 cbm
Weight of case	2.22 kg
EAN inner box	4038917124614
EAN case	4038917047975

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT