



Product specification

Screw cap micro tube, 0.5 ml, Biosphere® plus



Product information

Product number	72.730.217
----------------	------------

Product descriptions

Name	Screw cap micro tube, 0.5 ml, Biosphere® plus
Product number	72.730.217
Product description	Screw cap micro tube, working volume: 0.5 ml, skirted conical base, with knurling, transparent, cap: natural, cap attached, with printed writing space, with graduations, Biosphere® plus, 25 piece(s)/bag

Material & colours

Color of product	transparent
O-Ring material	Ethylene propylene diene rubber (EPDM)
Closure material	Polypropylene (PP)
Color of cap	natural

Size

Diameter	10.8 mm
Length of product	44 mm

Product characteristics

Base shape	skirted conical base
------------	----------------------

Label/Print	With print
Color of print/label	white
Graduation	yes
Knurls	yes
Cap	Retaining loop cap
Centrifugation max (RCF)	20000 x g

Purity & certification

Satisfies the requirement	IATA
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	250
Piece(s) / case	500
Piece(s) / pallet	49500
Height of inner box	141 mm
Width of inner box	191 mm
Depth of inner box	211 mm
Height of case	160 mm
Width of case	223 mm
Depth of case	398 mm
Volume of case	0.0142 cbm
Weight of case	1.39 kg
EAN inner box	4038917224017
EAN case	4038917224000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

