



Product specification

Screw cap micro tube, 0.5 ml, PCR Performance Tested



Product information

Product number	72.730.406
----------------	------------

Product descriptions

Name	Screw cap micro tube, 0.5 ml, PCR Performance Tested
Product number	72.730.406
Product description	Screw cap micro tube, working volume: 0.5 ml, skirted conical base, with knurling, transparent, cap: natural, cap assembled, with printed writing space, PCR Performance Tested, 100 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
O-Ring material	Ethylene propylene diene rubber (EPDM)
Closure material	Polypropylene (PP)
Color of cap	natural

Size

Diameter	10.8 mm
Length of product	44 mm

Product characteristics

Base shape	skirted conical base
Label/Print	With print
Color of print/label	white
Graduation	no
Knurls	yes
Cap	standard screw cap
Closure type	Screw cap
Centrifugation max (RCF)	20000 x g

Purity & certification

Satisfies the requirement	IATA
Product category	diagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order quantity	1000
Piece(s) / inner box	500
Piece(s) / case	1000
Piece(s) / pallet	110000
Height of inner box	141 mm
Width of inner box	191 mm
Depth of inner box	211 mm
Height of case	160 mm
Width of case	223 mm
Depth of case	398 mm
Volume of case	0.0142 cbm
Weight of case	2.04 kg
EAN inner box	4038917223980
EAN case	4038917223973

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support
and is valid without signature.



SARSTEDT