



Product specification

PCR single tube, 0.2 ml, Biosphere® plus, transparent, PP, flat cap



Product information

| | |
|----------------|--------|
| Product number | 72.737 |
|----------------|--------|

Product descriptions

| | |
|---------------------|--|
| Name | PCR single tube, 0.2 ml, Biosphere® plus, transparent, PP, flat cap |
| Product number | 72.737 |
| Product description | PCR single tube, 0.2 ml, Biosphere® plus, transparent, material: PP, with anti-contamination shield, flat cap, suitable for qPCR, 250 piece(s)/bag |

Material & colours

| | |
|------------------|--------------|
| Color of product | transparent |
| Profile | High Profile |

Size

| | |
|-------------------|----------|
| Diameter | 7.5 mm |
| Width of product | 7.5 mm |
| Length of product | 21.7 mm |
| Height of product | 20.65 mm |
| Max. volume | 200 µl |
| Sample volume | 0.2 ml |

Product characteristics

| | |
|----------------------------|-------------|
| Base shape | conical |
| Graduation | no |
| Cap | flat cap |
| Closure type | Snap-on lid |
| Autoclavable (max. °C) | 121 °C |
| Suitable for qPCR | yes |
| Centrifugation max (RCF) | 8000 x g |
| Centrifugation temperature | 22 °C |
| Centrifugation time | 5 Min |

Purity & certification

| | |
|------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD) |
| CE certified | CE - manufacturer's self-declaration |
| Purity standard | Biosphere plus |
| Sterilisation | Gas sterilisation |
| Pyrogen-free | yes |
| Batched | yes |

Packaging

| | |
|----------------------|---------------|
| Piece(s) / inner box | 1000 |
| Piece(s) / case | 2000 |
| Piece(s) / pallet | 288000 |
| Height of inner box | 97 mm |
| Width of inner box | 151 mm |
| Depth of inner box | 189 mm |
| Height of case | 145 mm |
| Width of case | 213 mm |
| Depth of case | 323 mm |
| Volume of case | 0.01 cbm |
| Weight of case | 0.72 kg |
| EAN inner box | 4038917083836 |
| EAN case | 4038917028370 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT