



Product specification

**PCR single tube, 0.2 ml, PCR Performance Tested,
transparent, PP, flat cap**



Product information

Product number	72.737.002
----------------	------------

Product descriptions

Name	PCR single tube, 0.2 ml, PCR Performance Tested, transparent, PP, flat cap
Product number	72.737.002
Product description	PCR single tube, 0.2 ml, PCR Performance Tested, transparent, material: PP, with anti-contamination shield, flat cap, suitable for qPCR, 500 piece(s)/minigrip bag

Material & colours

Color of product	transparent
Profile	High Profile

Size

Diameter	7.5 mm
Width of product	7.5 mm
Length of product	21.7 mm
Height of product	20.65 mm
Max. volume	200 µl
Sample volume	0.2 ml

Product characteristics

Base shape	conical
Graduation	no
Cap	flat cap
Closure type	Push cap
Autoclavable (max. °C)	121 °C
Suitable for qPCR	yes
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	648000
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm
Weight of case	0.53 kg
EAN inner box	4038917120104
EAN case	4038917017060

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



SARSTEDT

This document was prepared by EDP support and is valid without signature.