




Product specification

**PCR strip of 8, 100 µl, PCR Performance Tested,
transparent, PP, flat cap**



Product information

Product number	72.982.002
----------------	------------

Product descriptions

Name	PCR strip of 8, 100 µl, PCR Performance Tested, transparent, PP, flat cap
Product number	72.982.002
Product description	PCR strip of 8, up to 100 µl, Low Profile, PCR Performance Tested, transparent, material: PP, flat cap, cap enclosed, suitable for qPCR, 125 piece(s)/minigrip bag

Material & colours

Color of product	transparent
Profile	Low Profile

Size

Width of product	16.6 mm
Length of product	76.2 mm
Max. volume	100 µl
Sample volume	0.1 ml

Product characteristics

Base shape	conical
Graduation	no
Cap	flat cap
Closure type	Push cap
Autoclavable (max. °C)	121 °C
Suitable for qPCR	yes
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Piece(s) / inner box	125
Piece(s) / case	1250
Piece(s) / pallet	56250
Height of inner box	65 mm
Width of inner box	180 mm
Depth of inner box	220 mm
Height of case	198 mm
Width of case	300 mm
Depth of case	495 mm
Volume of case	0.0294 cbm
Weight of case	3.18 kg
EAN inner box	4038917395465
EAN case	4038917395458

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

