



## Product specification

**PCR strip of 8, 200 µl, PCR Performance Tested,  
transparent, PP**



### Product information

Product number	72.985.002
----------------	------------

### Product descriptions

Name	PCR strip of 8, 200 µl, PCR Performance Tested, transparent, PP
Product number	72.985.002
Product description	PCR strip of 8, up to 200 µl, High Profile, PCR Performance Tested, transparent, material: PP, without cap, suitable for qPCR, 120 piece(s)/minigrip bag

### Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Profile	High Profile

### Size

Diameter	5.4 mm
Width of product	8.1 mm
Length of product	70.9 mm
Height of product	20.7 mm
Max. volume	200 µl

Sample volume

0.2 ml

## Product characteristics

Base shape	conical
Graduation	no
Autoclavable (max. °C)	121 °C
Suitable for qPCR	yes
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min

## Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

## Packaging

Minimum order quantity	480
Piece(s) / inner box	120
Piece(s) / case	480
Piece(s) / pallet	155520
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm
Weight of case	0.6302 kg
EAN inner box	4038917124027
EAN case	4038917019002

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



# SARSTEDT

This document was prepared by EDP support and is valid without signature.