



Product specification

**PCR strip of 8, 200 µl, PCR Performance Tested,
transparent, PP, flat cap**



Product information

Product number	72.991.002
----------------	------------

Product descriptions

Name	PCR strip of 8, 200 µl, PCR Performance Tested, transparent, PP, flat cap
Product number	72.991.002
Product description	PCR strip of 8, up to 200 µl, High Profile, PCR Performance Tested, transparent, material: PP, with anti-contamination shield, flat cap, cap attached, suitable for qPCR, 120 piece(s)/minigrip bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Profile	High Profile

Size

Diameter	7.5 mm
Width of product	10.2 mm
Length of product	72.4 mm
Height of product	20.65 mm

Max. volume	200 µl
Sample volume	0.2 ml

Product characteristics

Base shape	conical
Graduation	no
Cap	flat cap
Closure type	Push cap
Autoclavable (max. °C)	121 °C
Suitable for qPCR	yes
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order quantity	480
Piece(s) / inner box	120
Piece(s) / case	480
Piece(s) / pallet	34560
Height of case	140 mm
Width of case	323 mm
Depth of case	358 mm
Volume of case	0.0162 cbm
Weight of case	0.9902 kg
EAN inner box	4038917124294
EAN case	4038917067546

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support
and is valid without signature.



SARSTEDT