



## Product specification

**Insert, suitable for tubes and S-Monovette® Ø 15 mm,  
transparent**



### Product information

Product number	73.1052
----------------	---------

### Product descriptions

Name	Insert, suitable for tubes and S-Monovette® Ø 15 mm, transparent
Product number	73.1052
Product description	Insert, suitable for tubes and S-Monovette® Ø 15 mm, material: PS, 1.5 ml, length: 23.5 mm, transparent, 500 piece(s)/bag

### Material & colours

Color of product	transparent
------------------	-------------

### Size

Length of product	23.5 mm
Diameter of opening	11 mm
Max. volume	1.5 ml

### Product characteristics

Base shape	skirted conical base
Graduation	no

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

## Packaging

Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	156000
Height of case	173 mm
Width of case	355 mm
Depth of case	390 mm
Volume of case	0.024 cbm
Weight of case	2.02 kg
EAN inner box	4038917085151
EAN case	4038917003834

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT