



Product specification

**Insert, suitable for tubes and S-Monovette® Ø 15 mm,
transparent**



Product information

Product number	73.1052
----------------	---------

Product descriptions

Name	Insert, suitable for tubes and S-Monovette® Ø 15 mm, transparent
Product number	73.1052
Product description	Insert, suitable for tubes and S-Monovette® Ø 15 mm, material: PS, length: 23.5 mm, transparent, 500 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

Size

Length of product	23.5 mm
Diameter of opening	11 mm

Product characteristics

Base shape	skirted conical base
Graduation	no

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order quantity	2000
Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	156000
Height of case	130 mm
Width of case	323 mm
Depth of case	428 mm
Volume of case	0.018 cbm
Weight of case	2.02 kg
EAN inner box	4038917085151
EAN case	4038917003834

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT