



## Product specification

**Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, PP, transparent**



## Product information

Product number	75.562.005
----------------	------------

## Product descriptions

Name	Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, PP, transparent
Product number	75.562.005
Product description	Container with screw cap, collection and storage of urine, max. working volume: 100 ml, (ØxH): 57 x 76 mm, Ø opening: 62 mm, transparent, graduated, material: PP, screw cap, cap: white, cap assembled, 5 piece(s)/bag

## Material & colours

Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	white

## Size

Diameter	57 mm
Length excluding cap	72 mm
Height of product	76 mm
Diameter of opening	62 mm

Max. volume	100 ml
Sample volume	100 ml

## Product characteristics

Graduation	yes
Closure type	Screw cap

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

## Packaging

Piece(s) / inner box	200
Piece(s) / case	200
Piece(s) / pallet	5600
Height of case	255 mm
Width of case	400 mm
Depth of case	590 mm
Volume of case	0.0602 cbm
Weight of case	3.74 kg
EAN inner box	4038917113595
EAN case	4038917014021

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT