



Product specification

Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, with light protection, PP, brown



Product information

Product number	75.562.011
----------------	------------

Product descriptions

Name	Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, with light protection, PP, brown
Product number	75.562.011
Product description	Container with screw cap, collection and storage of urine, max. working volume: 100 ml, (ØxH): 57 x 76 mm, Ø opening: 62 mm, brown, with light protection, graduated, material: PP, screw cap, cap: brown, cap enclosed, 50 piece(s)/bag

Material & colours

Color of product	brown
Closure material	High Density Polyethylene (HD-PE)
Color of cap	brown

Size

Diameter	57 mm
Length excluding cap	72 mm
Height of product	76 mm

Diameter of opening	62 mm
Max. volume	100 ml
Sample volume	100 ml

Product characteristics

Graduation	yes
Light protection	yes
Closure type	Screw cap

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	250
Piece(s) / pallet	7500
Height of case	330 mm
Width of case	323 mm
Depth of case	443 mm
Volume of case	0.0472 cbm
Weight of case	4.54 kg
EAN inner box	4038917113618
EAN case	4038917014038

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT