



Product specification

**Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, PP,
with safety label, transparent**



Product information

Product number	75.562.105
----------------	------------

Product descriptions

Name	Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, PP, with safety label, transparent
Product number	75.562.105
Product description	Container with screw cap, collection and storage of urine, max. working volume: 100 ml, (ØxH): 57 x 76 mm, Ø opening: 62 mm, transparent, graduated, material: PP, screw cap, cap: white, with safety label, cap assembled, sterile, 5 piece(s)/bag

Material & colours

Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	white

Size

Diameter	57 mm
Length excluding cap	72 mm
Height of product	76 mm

Diameter of opening	62 mm
Max. volume	100 ml
Sample volume	100 ml

Product characteristics

Label/Print	With safety label
Color of print/label	white/red
Graduation	yes
Closure type	Screw cap

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Piece(s) / case	200
Piece(s) / pallet	5600
Height of case	255 mm
Width of case	400 mm
Depth of case	590 mm
Volume of case	0.0602 cbm
Weight of case	3.78 kg
EAN inner box	4038917113632
EAN case	4038917014052

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT