



Product specification

Multi-purpose container, max. volume: 250 ml, (LxØ): 78 x 70 mm, graduated, PP, transparent



Product information

Product number	75.9922.534
----------------	-------------

Product descriptions

Name	Multi-purpose container, max. volume: 250 ml, (LxØ): 78 x 70 mm, graduated, PP, transparent
Product number	75.9922.534
Product description	Multi-purpose container, max. volume: 250 ml, (LxØ): 78 x 70 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, sterile, 80 piece(s)/bag

Material & colours

Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	yellow

Size

Diameter	70 mm
Length of product	78 mm
Height of product	78 mm
Max. volume	250 ml
Sample volume	250 ml

Product characteristics

Base shape	Flat base
Graduation	yes
Closure type	Screw cap

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Piece(s) / inner box	80
Piece(s) / case	240
Piece(s) / pallet	3840
Height of inner box	498 mm
Width of inner box	450 mm
Depth of inner box	646 mm
Height of case	449 mm
Width of case	498 mm
Depth of case	646 mm
Volume of case	0.1444 cbm
Weight of case	10.865 kg
EAN inner box	4066916011215
EAN case	4066916011222

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT