



## Product specification

**Multi-purpose container, 250 ml, (LxØ): 78 x 70 mm,  
graduated, PP, transparent**



## Product information

Product number	75.9922.534
----------------	-------------

## Product descriptions

Name	Multi-purpose container, 250 ml, (LxØ): 78 x 70 mm, graduated, PP, transparent
Product number	75.9922.534
Product description	Multi-purpose container, max. working volume: 250 ml, (LxØ): 78 x 70 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, sterile, 240 piece(s)/bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	yellow

## Size

Diameter	70 mm
Length of product	78 mm
Height of product	78 mm
Sample volume	250 ml

## Product characteristics

Base shape	Flat base
Graduation	yes
Closure type	Screw cap

## Purity & certification

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

## Packaging

Minimum order quantity	240
Piece(s) / inner box	240
Piece(s) / case	240
Piece(s) / pallet	1440
Height of case	498 mm
Width of case	450 mm
Depth of case	646 mm
Volume of case	0.1448 cbm
Weight of case	10.9301 kg
EAN inner box	4038917152129
EAN case	4038917350907

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT