



Product specification

Multi-purpose container, max. volume: 70 ml, (LxØ): 55 x 44 mm, graduated, PP, transparent, with pap



Product information

Product number	75.9922.721
----------------	-------------

Product descriptions

Name	Multi-purpose container, max. volume: 70 ml, (LxØ): 55 x 44 mm, graduated, PP, transparent, with paper label
Product number	75.9922.721
Product description	Multi-purpose container, max. volume: 70 ml, (LxØ): 55 x 44 mm, Ø opening: 44 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, with paper label, sterile, 1 piece(s)/blister

Material & colours

Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	yellow

Size

Diameter	44 mm
Length excluding cap	55 mm
Length of product	55 mm
Height of product	55 mm

Diameter of opening	44 mm
Max. volume	70 ml

Product characteristics

Base shape	Flat base
Label/Print	With paper label
Graduation	yes
Closure type	Screw cap

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Piece(s) / inner box	200
Piece(s) / case	200
Piece(s) / pallet	3200
Height of case	405 mm
Width of case	375 mm
Depth of case	548 mm
Volume of case	0.0832 cbm
Weight of case	3.5 kg
EAN inner box	4038917152525
EAN case	4038917347181

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT