



## Product specification

**Multi-purpose container, max. volume: 70 ml, (LxØ): 55 x 44 mm, graduated, PP, transparent**



## Product information

Product number	75.9922.744
----------------	-------------

## Product descriptions

Name	Multi-purpose container, max. volume: 70 ml, (LxØ): 55 x 44 mm, graduated, PP, transparent
Product number	75.9922.744
Product description	Multi-purpose container, max. volume: 70 ml, (LxØ): 55 x 44 mm, Ø opening: 44 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, 250 piece(s)/bag

## Material & colours

Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	yellow

## Size

Diameter	44 mm
Length excluding cap	55 mm
Length of product	55 mm
Height of product	55 mm
Diameter of opening	44 mm

Max. volume

70 ml

## Product characteristics

Base shape	Flat base
Graduation	yes
Closure type	Screw cap

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE
Batched	yes

## Packaging

Piece(s) / inner box	250
Piece(s) / case	500
Piece(s) / pallet	8000
Height of case	373 mm
Width of case	408 mm
Depth of case	548 mm
Volume of case	0.0834 cbm
Weight of case	7 kg
EAN inner box	4038917152754
EAN case	4038917027793

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT