




## Product specification

**Urine container, 70 ml, (LxØ): 55 x 44 mm, graduated, PP, transparent, with paper label**



## Product information

Product number	75.9922.745
----------------	-------------

## Product descriptions

Name	Urine container, 70 ml, (LxØ): 55 x 44 mm, graduated, PP, transparent, with paper label
Product number	75.9922.745
Product description	Urine container, max. working volume: 70 ml, (LxØ): 55 x 44 mm, Ø opening: 44 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, with paper label, label/print: green, 250 piece(s)/bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	yellow

## Size

Diameter	44 mm
Length excluding cap	55 mm
Length of product	55 mm

Height of product	55 mm
Diameter of opening	44 mm
Max. volume	70 ml
Sample volume	70 ml

## Product characteristics

Base shape	Flat base
Label/Print	With paper label
Color of print/label	green
Graduation	yes
Closure type	Screw cap

## Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

## Packaging

Minimum order quantity	500
Piece(s) / inner box	250
Piece(s) / case	500
Piece(s) / pallet	8000
Height of case	373 mm
Width of case	408 mm
Depth of case	548 mm
Volume of case	0.0834 cbm
Weight of case	6.6 kg
EAN inner box	4038917152761
EAN case	4038917027809

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



# SARSTEDT