



## Product specification

**Multi-purpose container, 420 ml, (LxØ): 150 x 70 mm,  
graduated, PP, transparent**



## Product information

Product number	75.9922.813
----------------	-------------

## Product descriptions

Name	Multi-purpose container, 420 ml, (LxØ): 150 x 70 mm, graduated, PP, transparent
Product number	75.9922.813
Product description	Multi-purpose container, max. working volume: 420 ml, (LxØ): 150 x 70 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, sterile, 40 piece(s)/bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Tube material	Polypropylene (PP)
Closure material	High Density Polyethylene (HD-PE)
Color of cap	yellow

## Size

Diameter	70 mm
Length of product	150 mm
Height of product	150 mm

Sample volume

420 ml

## Product characteristics

Base shape	Flat base
Graduation	yes
Cap	standard screw cap
Closure type	Screw cap

## Purity & certification

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Latex-free	yes
Batched	yes

## Packaging

Minimum order quantity	80
Piece(s) / inner box	40
Piece(s) / case	80
Piece(s) / pallet	1280
Height of case	408 mm
Width of case	373 mm
Depth of case	548 mm
Volume of case	0.0834 cbm
Weight of case	6.02 kg
EAN inner box	4038917153331
EAN case	4038917350938

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT