



Product specification

Snap-on lid, without spout, PE, red, 1,000 piece(s)/bag



Product information

Product number	76.1414
----------------	---------

Product descriptions

Name	Snap-on lid, without spout, PE, red, 1,000 piece(s)/Bag
Product number	76.1414
Product description	Snap-on lid, without spout, suitable for 125-ml container 75.1355, material: PE, red, 1,000 piece(s)/bag

Material & colours

Color of product	red
------------------	-----

Product characteristics

Design	without spout
Closure type	Snap-on lid

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE supplier
Batched	yes

Packaging

Piece(s) / inner box	1000
----------------------	------

Piece(s) / case	1000
Piece(s) / pallet	24000
Height of inner box	195 mm
Width of inner box	190 mm
Depth of inner box	380 mm
Height of case	285 mm
Width of case	390 mm
Depth of case	590 mm
Volume of case	0.0656 cbm
Weight of case	3.14 kg
EAN inner box	4038917086554
EAN case	4042836001499

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT