



Product specification

Absorbent liner, suitable for Mailing container 126 x 30 mm, (LxW): 75 x 90 mm



Product information

Product number	78.670
----------------	--------

Product descriptions

Name	Absorbent liner, suitable for mailing container 126 x 30 mm, (LxW): 75 x 90 mm
Product number	78.670
Product description	Absorbent liner, suitable for mailing container 126 x 30 mm, (LxW): 75 x 90 mm, 1,000 piece(s)/case

Material & colours

Product material	Cellulose
------------------	-----------

Size

Width of product	90 mm
Length of product	75 mm

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order quantity	4000
Piece(s) / inner box	1000
Piece(s) / case	4000
Piece(s) / pallet	72000
Height of case	360 mm
Width of case	360 mm
Depth of case	515 mm
Volume of case	0.0667 cbm
Weight of case	4.36 kg
EAN inner box	4038917377492
EAN case	4038917159937

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT