



## Product specification

**Faeces tube, with spoon, screw cap, (LxØ): 107 x 25 mm, transparent, sterile**



### Product information

Product number	80.622
----------------	--------

### Product descriptions

Name	Faeces tube, with Spoon, Screw cap, (LxØ): 107 x 25 mm, transparent, sterile
Product number	80.622
Product description	Faeces tube, with spoon, transparent, screw cap, Screw cap with stool spoon holder, cap: brown, (LxØ) without cap: 107 x 25 mm, sterile, 50 piece(s)/bag

### Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)
Stem material	Polystyrene (PS)
Closure material	High Density Polyethylene (HD-PE)
Color of cap	brown

### Size

Diameter	25 mm
Length excluding cap	107 mm
Length of product	111 mm

## Product characteristics

Design	Spoon
Base shape	Flat-conical base
Cap	Screw cap with stool spoon holder
Closure type	Screw cap

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Piece(s) / inner box	50
Piece(s) / case	250
Piece(s) / pallet	10000
Height of inner box	50 mm
Width of inner box	395 mm
Depth of inner box	435 mm
Height of case	165 mm
Width of case	398 mm
Depth of case	598 mm
Volume of case	0.0393 cbm
Weight of case	4.29 kg
EAN inner box	4038917082976
EAN case	4066916010805

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT