



## Product specification

**Faeces tube, with spoon, screw cap, (LxØ): 107 x 25 mm, transparent**



### Product information

Product number	80.622.111
----------------	------------

### Product descriptions

Name	Faeces tube, with spoon, screw cap, (LxØ): 107 x 25 mm, transparent
Product number	80.622.111
Product description	Faeces tube, with spoon, transparent, screw cap, Screw cap with stool spoon, enclosed, cap: brown, (LxØ) without cap: 107 x 25 mm, 250 piece(s)/bag

### Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)
Stem material	Polystyrene (PS)
Closure material	High Density Polyethylene (HD-PE)
Color of cap	brown

### Size

Diameter	25 mm
Length excluding cap	107 mm
Length of product	111 mm

## Product characteristics

Design	Spoon
Base shape	Flat-conical base
Cap	Screw cap with stool spoon
Closure type	Screw cap

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

## Packaging

Piece(s) / inner box	250
Piece(s) / case	500
Piece(s) / pallet	12000
Height of case	415 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0657 cbm
Weight of case	7.14 kg
EAN inner box	4038917115346
EAN case	4038917014533

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT