



Product specification

Faeces tube, with spoon, screw cap, (LxØ): 101 x 16.5 mm, transparent, sterile



Product information

Product number	80.623
----------------	--------

Product descriptions

Name	Faeces tube, with spoon, screw cap, (LxØ): 101 x 16.5 mm, transparent, sterile
Product number	80.623
Product description	Faeces tube, with spoon, transparent, screw cap, Screw cap with stool spoon, assembled, cap: brown, (LxØ) without cap: 101 x 16.5 mm, outer diameter of cap: 19.5 mm, sterile, 500 piece(s)/bag

Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)
Stem material	Polystyrene (PS)
Closure material	High Density Polyethylene (HD-PE)
Color of cap	brown

Size

Diameter	16.5 mm
Length excluding cap	101 mm

Length of product	107 mm
-------------------	--------

Product characteristics

Design	Spoon
Base shape	Round base
Cap	Screw cap with stool spoon
Closure type	Screw cap

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	500
Piece(s) / pallet	22000
Height of case	155 mm
Width of case	348 mm
Depth of case	568 mm
Volume of case	0.0306 cbm
Weight of case	3.89 kg
EAN inner box	4038917082983
EAN case	4038917001960

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT