



## Product specification

**Faeces tube, with spoon, screw cap, (LxØ): 101 x 16.5 mm, transparent, sterile**



### Product information

|                |        |
|----------------|--------|
| Product number | 80.623 |
|----------------|--------|

### Product descriptions

|                     |   |
|---------------------|---|
| Name                | Faeces tube, with spoon, screw cap, (LxØ): 101 x 16.5 mm, transparent, sterile  |
| Product number      | 80.623  |
| Product description | Faeces tube, with spoon, transparent, screw cap, Screw cap with stool spoon, assembled, cap: brown, (LxØ) without cap: 101 x 16.5 mm, sterile, 500 piece(s)/bag |

### Material & colours

|                  |                                   |
|------------------|-----------------------------------|
| Color of product | transparent                       |
| Tube material    | Polypropylene (PP)                |
| Stem material    | Polystyrene (PS)                  |
| Closure material | High Density Polyethylene (HD-PE) |
| Color of cap     | brown                             |

### Size

|                      |         |
|----------------------|---------|
| Diameter             | 16.5 mm |
| Length excluding cap | 101 mm  |
| Length of product    | 107 mm  |

## Product characteristics

|              |                            |
|--------------|----------------------------|
| Design       | Spoon                      |
| Base shape   | Round base                 |
| Cap          | Screw cap with stool spoon |
| Closure type | Screw cap                  |

## Purity & certification

|                           |                       |
|---------------------------|-----------------------|
| Satisfies the requirement | IATA, ADR             |
| Product category          | diagnostic device, CE |
| CE certified              | CE                    |
| Purity standard           | sterile               |
| Sterilisation             | Electron irradiation  |
| Batched                   | yes                   |

## Packaging

|                        |               |
|------------------------|---------------|
| Minimum order quantity | 500           |
| Piece(s) / inner box   | 500           |
| Piece(s) / case        | 500           |
| Piece(s) / pallet      | 22000         |
| Height of case         | 155 mm        |
| Width of case          | 348 mm        |
| Depth of case          | 568 mm        |
| Volume of case         | 0.0306 cbm    |
| Weight of case         | 3.89 kg       |
| EAN inner box          | 4038917082983 |
| EAN case               | 4038917001960 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT