



## Product specification

**Faeces tube, with special faeces spoon (1 ml), screw cap,  
(LxØ): 101 x 16.5 mm, transparent**



## Product information

Product number	80.623.022
----------------	------------

## Product descriptions

Name	Faeces tube, with special faeces spoon (1 ml), screw cap, (LxØ): 101 x 16.5 mm, transparent
Product number	80.623.022
Product description	Faeces tube, with special faeces spoon (1 ml), With enclosed stirring rod, transparent, screw cap, cap: brown, (LxØ) without cap: 101 x 16.5 mm, with paper label, 250 piece(s)/bag

## Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)
Closure material	High Density Polyethylene (HD-PE)
Color of cap	brown

## Size

Diameter	16.5 mm
Length excluding cap	101 mm
Length of product	107 mm

## Product characteristics

Design	special faeces spoon (1 ml)
Label/Print	With paper label
Color of print/label	black
Closure type	Screw cap

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

## Packaging

Piece(s) / inner box	250
Piece(s) / case	1000
Piece(s) / pallet	24000
Height of case	415 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0657 cbm
Weight of case	8.62 kg
EAN inner box	4038917115407
EAN case	4038917014601

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT