



Product specification

Forensic swab, round, in the tube with ventilation membrane, ISO 18385, 95 mm, cotton



Product information

Product number	80.629.001
----------------	------------

Product descriptions

Name	Forensic swab, round, in the tube with ventilation membrane, ISO 18385, 95 mm, cotton
------	---

Product number	80.629.001
----------------	------------

Product description	Forensic swab, round, in the tube with ventilation membrane, with paper label, ETO treated, ISO 18385, passive drying, stem length: 95 mm, stem material: wood, swab tip round, swab material: cotton, tube: PP, blister packaging aluminium free, 1 piece(s)/blister
---------------------	---

Material & colours

Color of product	multi-coloured
------------------	----------------

Tube material	Polypropylene (PP)
---------------	--------------------

Stem material	Wood
---------------	------

Swab material	Cotton
---------------	--------

Closure material	High Density Polyethylene (HD-PE)
------------------	-----------------------------------

Color of cap	white
--------------	-------

Size

Diameter	15.35 mm
Length excluding cap	106 mm
Length including cap	97.6 µm
Width of product	18.05 mm
Length of product	107 mm
Height of product	106.8 mm
Length including cap and plunger	106.8 µm
Length of the stem	95 mm

Product characteristics

Design	in the tube with ventilation membrane
Label/Print	With paper label
Cap	standard screw cap
Closure type	Screw cap

Purity & certification

Satisfies the requirement	ISO 18385
Product category	no medical device IVD
Purity standard	ETO treated
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order quantity	800
Type of inner packaging	Box
Piece(s) / inner box	100
Piece(s) / case	800
Piece(s) / pallet	20000
Height of inner box	171 mm
Width of inner box	151 mm
Depth of inner box	211 mm
Height of case	335 mm
Width of case	363 mm
Depth of case	437 mm
Volume of case	0.0531 cbm
Weight of case	7.2904 kg
EAN inner box	4038917276887
EAN case	4038917428941

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT