



Product specification

Faeces tube, with spoon, screw cap, (LxØ): 76 x 20 mm, transparent, sterile



Product information

Product number	80.734.001
----------------	------------

Product descriptions

Name	Faeces tube, with spoon, screw cap, (LxØ): 76 x 20 mm, transparent, sterile
Product number	80.734.001
Product description	Faeces tube, with spoon, transparent, screw cap, Screw cap with stool spoon, assembled, cap: brown, (LxØ) without cap: 76 x 20 mm, sterile, 100 piece(s)/bag

Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)
Stem material	Polystyrene (PS)
Closure material	High Density Polyethylene (HD-PE)
Color of cap	brown

Size

Diameter	20 mm
Length excluding cap	76 mm
Length of product	77 mm

Product characteristics

Design	Spoon
Base shape	Flat-conical base
Cap	Screw cap with stool spoon
Closure type	Screw cap

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	20000
Height of case	165 mm
Width of case	398 mm
Depth of case	598 mm
Volume of case	0.0393 cbm
Weight of case	4.06 kg
EAN inner box	4038917119801
EAN case	4038917016810

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT