




Product specification

Faeces tube, screw cap, (LxØ): 76 x 20 mm, transparent, sterile



Product information

Product number	80.734.311
----------------	------------

Product descriptions

Name	Faeces tube, Screw cap, (LxØ): 76 x 20 mm, transparent, sterile
Product number	80.734.311
Product description	Faeces tube, transparent, screw cap, cap: brown, (LxØ) without cap: 76 x 20 mm, with paper label, sterile, 100 piece(s)/bag

Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)
Closure material	High Density Polyethylene (HD-PE)
Color of cap	brown

Size

Diameter	20 mm
Length excluding cap	76 mm
Length of product	77 mm

Product characteristics

Label/Print	With paper label
Color of print/label	black
Closure type	Screw cap

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	20000
Height of case	165 mm
Width of case	398 mm
Depth of case	598 mm
Volume of case	0.0393 cbm
Weight of case	4.49 kg
EAN inner box	4038917119962
EAN case	4038917016964

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT