



Product specification

Petri dish, 92 x 16 mm, transparent, with ventilation cams



Product information

Product number	82.1473.001
----------------	-------------

Product descriptions

Name	Petri dish, 92 x 16 mm, transparent, with ventilation cams
------	--

Product number	82.1473.001
----------------	-------------

Product description	Petri dish, (ØxH): 92 x 16 mm, material: PS, transparent, with ventilation cams, sterile, 20 piece(s)/tube bag
---------------------	--

Material & colours

Product material	Polystyrene (PS)
------------------	------------------

Color of product	transparent
------------------	-------------

Size

Diameter	92 mm
----------	-------

Height of product	16 mm
-------------------	-------

Max. volume	80 ml
-------------	-------

Product characteristics

Ventilation cams	yes
------------------	-----

Minimum temperature	0 °C
---------------------	------

Maximum temperature	80 °C
---------------------	-------

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Latex-free	yes
Batched	yes
Safety data sheet	no

Packaging

Minimum order quantity	480
Type of inner packaging	Tube bag
Piece(s) / inner box	480
Piece(s) / case	480
Piece(s) / pallet	9600
Height of case	360 mm
Width of case	387 mm
Depth of case	584 mm
Volume of case	0.0814 cbm
Weight of case	8.7902 kg
EAN inner box	4038917129329
EAN case	4038917022125

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT