



Product specification

Petri dish, 92 x 16 mm, yellow, with ventilation cams



Product information

Product number	82.1473.040
----------------	-------------

Product descriptions

Name	Petri dish, 92 x 16 mm, yellow, with ventilation cams
------	---

Product number	82.1473.040
----------------	-------------

Product description	Petri dish, (ØxH): 92 x 16 mm, material: PS, yellow, with ventilation cams, 20 piece(s)/tube bag
---------------------	--

Material & colours

Color of product	yellow
------------------	--------

Size

Diameter	92 mm
----------	-------

Height of product	16 mm
-------------------	-------

Max. volume	80 ml
-------------	-------

Product characteristics

Ventilation cams	yes
------------------	-----

Minimum temperature	0 °C
---------------------	------

Maximum temperature	80 °C
---------------------	-------

Purity & certification

Product category	not regulated
------------------	---------------

CE certified	no CE mark
Latex-free	yes
Batched	yes
Safety data sheet	no

Packaging

Type of inner packaging	Tube bag
Piece(s) / inner box	480
Piece(s) / case	480
Piece(s) / pallet	9600
Height of case	350 mm
Width of case	387 mm
Depth of case	584 mm
Volume of case	0.0791 cbm
Weight of case	9.2501 kg
EAN inner box	4038917129374
EAN case	4038917022156

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT