



## Product specification

**Micro test plate, 96 well, slip-on lid, flat base, PS,  
transparent**



## Product information

Product number	82.1581.001
----------------	-------------

## Product descriptions

Name	Micro test plate, 96 Well, Slip-on lid, Flat base, PS, transparent
Product number	82.1581.001
Product description	Micro test plate, 96 well, slip-on lid, flat base, lid assembled, material: PS, transparent, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 1 piece(s)/blister

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Width of product	85.2 mm
Length of product	127.8 mm
Height of product	16.55 mm
Max. volume	390 µl

## Product characteristics

Number of wells	96
-----------------	----

Base shape	Flat base
Label/Print	With plastic label
Vessel type	Plate
Closure type	Slip-on lid

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation	Gamma irradiation
Batched	yes

## Packaging

Piece(s) / case	50
Piece(s) / pallet	3200
Height of case	205 mm
Width of case	297 mm
Depth of case	332 mm
Volume of case	0.0202 cbm
Weight of case	3.54 kg
EAN inner box	4038917131056
EAN case	4038917347952

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT