



Product specification

ELISA plate, 96 well, flat base, PS, black, Medium Binding



Product information

Product number	82.1581.120
----------------	-------------

Product descriptions

Name	ELISA plate, 96 well, flat base, PS, black, Medium Binding
Product number	82.1581.120
Product description	ELISA plate, 96 well, flat base, without lid, material: PS, black, Medium Binding, free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic, 25 piece(s)/bag

Material & colours

Color of product	black
------------------	-------

Size

Width of product	85.2 mm
Length of product	127.8 mm
Height of product	14.35 mm
Max. volume	390 µl

Product characteristics

Number of wells	96
Base shape	Flat base
Label/Print	With print

Vessel type

Plate

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	25
Piece(s) / case	50
Piece(s) / pallet	2800
Height of case	365 mm
Width of case	232 mm
Depth of case	279 mm
Volume of case	0.0236 cbm
Weight of case	2.86 kg
EAN inner box	4038917365338
EAN case	4038917365321

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT