



## Product specification

# ELISA plate, 96 well, round base, PS, transparent, High Binding



## Product information

Product number	82.1582.200
----------------	-------------

## Product descriptions

Name	ELISA plate, 96 well, round base, PS, transparent, High Binding
Product number	82.1582.200
Product description	ELISA plate, 96 well, round base, without lid, material: PS, transparent, High Binding, free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic, 25 piece(s)/bag

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Width of product	85.2 mm
Length of product	127.8 mm
Height of product	14.35 mm
Max. volume	310 µl

## Product characteristics

Number of wells	96
-----------------	----

Base shape	Round base
Label/Print	With print
Vessel type	Plate

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic
Pyrogen-free	yes
Batched	yes

## Packaging

Piece(s) / inner box	25
Piece(s) / case	50
Piece(s) / pallet	2000
Height of case	365 mm
Width of case	232 mm
Depth of case	279 mm
Volume of case	0.0236 cbm
Weight of case	2.54 kg
EAN inner box	4038917356220
EAN case	4038917356213

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT