



Product specification

Syringe filter, Filtropur S, PES, pore size: 0.45 µm, for clear filtration



Product information

Product number	83.1826
----------------	---------

Product descriptions

Name	Syringe filter, Filtropur S, PES, pore size: 0.45 µm, for clear filtration
Product number	83.1826
Product description	Syringe filter, Filtropur S, membrane: PES, filtration surface: 6.2 cm ² , pore size: 0.45 µm, for clear filtration, connection: Luer lock, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Material & colours

Product material	Methacrylate-butadiene-styrene (MBS)
Color of product	yellow
Membrane material	Polyethersulfone (PES)

Size

Diameter	33 mm
----------	-------

Product characteristics

Filter	yes
Area of the filter	6.2 cm ²

Pore size	0.45 µm
Connection	Luer lock
Example of use	for clear filtration

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	50
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	12000
Height of case	530 mm
Width of case	285 mm
Depth of case	375 mm
Volume of case	0.0566 cbm
Weight of case	3.1 kg
EAN inner box	4038917087742
EAN case	4038917006026

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT