



Product specification

Syringe filter, Filtropur S, PES, pore size: 0.45 µm, for clear filtration



Product information

| | |
|----------------|---------|
| Product number | 83.1826 |
|----------------|---------|

Product descriptions

| | |
|---------------------|--|
| Name | Syringe filter, Filtropur S, PES, pore size: 0.45 µm, for clear filtration |
| Product number | 83.1826 |
| Product description | Syringe filter, Filtropur S, membrane: PES, filtration surface: 6.2 cm², pore size: 0.45 µm, for clear filtration, connection: Luer lock, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister |

Material & colours

| | |
|-------------------|------------------------|
| Color of product | yellow |
| Membrane material | Polyethersulfone (PES) |

Size

| | |
|----------|-------|
| Diameter | 33 mm |
|----------|-------|

Product characteristics

| | |
|--------------------|---------|
| Filter | yes |
| Area of the filter | 6.2 cm² |
| Pore size | 0.45 µm |

| | |
|----------------|----------------------|
| Connection | Luer lock |
| Example of use | for clear filtration |

Purity & certification

| | |
|------------------|--|
| Product category | not regulated |
| CE certified | no CE mark |
| Purity standard | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation | Gamma irradiation |
| Pyrogen-free | yes |
| Batched | yes |

Packaging

| | |
|----------------------|---------------|
| Piece(s) / inner box | 50 |
| Piece(s) / case | 500 |
| Piece(s) / pallet | 12000 |
| Height of inner box | 76 mm |
| Width of inner box | 140 mm |
| Depth of inner box | 140 mm |
| Height of case | 280 mm |
| Width of case | 365 mm |
| Depth of case | 525 mm |
| Volume of case | 0.0537 cbm |
| Weight of case | 3.1 kg |
| EAN inner box | 4038917087742 |
| EAN case | 4038917006026 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT