



Product specification

Coverslips, Ø: 25 mm



Product information

Product number	83.1840
----------------	---------

Product descriptions

Name	Coverslips, Ø: 25 mm
Product number	83.1840
Product description	Coverslips, for adherent cells, Ø: 25 mm, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 25 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	25 mm
----------	-------

Product characteristics

Type of surface treatment	TC-treated
Cell type	adherent cells

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation

Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	200
Piece(s) / case	2000
Piece(s) / pallet	72000
Height of inner box	60 mm
Width of inner box	110 mm
Depth of inner box	145 mm
Height of case	144 mm
Width of case	310 mm
Depth of case	335 mm
Volume of case	0.015 cbm
Weight of case	1.8 kg
EAN inner box	4038917087841
EAN case	4038917006095

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT