

Product specification

Coverslips, Ø: 22 mm



Product information

| | |
|--------------|-------------|
| Order number | 83.1840.001 |
|--------------|-------------|

Product descriptions

| | |
|---------------------|--|
| Name | Coverslips, Ø: 22 mm |
| Product number | 83.1840.001 |
| Product description | Coverslips, cell type: adherent, Ø: 22 mm, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 25 piece(s)/bag |

Product characteristics

| | |
|---------------------------|---------------|
| Type of surface treatment | TC-treated |
| Suitable for | 6-well plates |
| Cell type | Adherent |

Size

| | |
|----------|-------|
| Diameter | 22 mm |
|----------|-------|

Material & colours

| | |
|------------------|---|
| Product material | Polyethylene terephthalate, Glycol-modified (PET-G) |
| Color of product | transparent |

Purity & certification

| | |
|------------------|--|
| Product category | no medical device IVD |
| Purity standard | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation | Gamma irradiation |
| Pyrogen-free | yes |
| Batched | yes |

Packaging

| | |
|----------------------|---------------|
| Piece(s) / inner box | 200 |
| Piece(s) / case | 2000 |
| Piece(s) / pallet | 72000 |
| Volume of case | 0.0175 cbm |
| Weight of case | 1.8 kg |
| EAN inner box | 4038917134040 |
| EAN case | 4038917025386 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

