




## Product specification

**Tissue culture dish, (ØxH): 35 x 10 mm, surface: Standard**



## Product information

Product number	83.3900
----------------	---------

## Product descriptions

Name	Tissue culture dish, (ØxH): 35 x 10 mm, Surface: Standard
Product number	83.3900
Product description	Tissue culture dish, (ØxH): 35 x 10 mm, material: PS, surface: Standard, for adherent cells, coding colour: red, TC Tested, 10 piece(s)/bag

## Material & colours

Closure material	Polystyrene (PS)
------------------	------------------

## Size

Diameter	35 mm
Height of product	10 mm

## Product characteristics

Growth area	9.4 cm <sup>2</sup>
Base shape	SUREGrip Design
Vessel type	Dish
Closure type	Slip-on lid
Cell type	adherent cells

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Piece(s) / inner box	10
Piece(s) / case	500
Piece(s) / pallet	30000
Height of case	255 mm
Width of case	257 mm
Depth of case	322 mm
Volume of case	0.0211 cbm
Weight of case	2.59 kg
EAN inner box	4038917337953
EAN case	4038917365789

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT