




Product specification

Tissue culture dish, (ØxH): 35 x 10 mm, surface: Cell+



Product information

Product number	83.3900.300
----------------	-------------

Product descriptions

Name	Tissue culture dish, (ØxH): 35 x 10 mm, Surface: Cell+
Product number	83.3900.300
Product description	Tissue culture dish, (ØxH): 35 x 10 mm, material: PS, surface: Cell+, for challenging adherent cells, coding colour: yellow, TC Tested, 10 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Closure material	Polystyrene (PS)

Size

Diameter	35 mm
Height of product	10 mm

Product characteristics

Growth area	9.4 cm ²
Base shape	SUREGrip Design
Vessel type	Dish
Closure type	Slip-on lid

Cell type	challenging adherent cells
-----------	----------------------------

Purity & certification

Product category	no medical device IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	10
Piece(s) / case	500
Piece(s) / pallet	30000
Height of case	245 mm
Width of case	257 mm
Depth of case	322 mm
Volume of case	0.0203 cbm
Weight of case	2.56 kg
EAN inner box	4038917338042
EAN case	4038917365802

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT