




Product specification

Tissue culture dish, (ØxH): 60 x 15 mm, surface: Standard



Product information

Product number	83.3901
----------------	---------

Product descriptions

Name	Tissue culture dish, (ØxH): 60 x 15 mm, Surface: Standard
Product number	83.3901
Product description	Tissue culture dish, (ØxH): 60 x 15 mm, material: PS, surface: Standard, for adherent cells, coding colour: red, TC Tested, 10 piece(s)/bag

Material & colours

Closure material	Polystyrene (PS)
------------------	------------------

Size

Diameter	60 mm
Height of product	15 mm

Product characteristics

Growth area	22.1 cm ²
Base shape	SUREGrip Design
Vessel type	Dish
Closure type	Slip-on lid
Cell type	adherent cells

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	10
Piece(s) / case	500
Piece(s) / pallet	12000
Height of case	380 mm
Width of case	332 mm
Depth of case	332 mm
Volume of case	0.0419 cbm
Weight of case	5.02 kg
EAN inner box	4038917338073
EAN case	4038917338066

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT