




## Product specification

**Tissue culture dish, (ØxH): 60 x 15 mm, surface: Standard**



## Product information

Product number	83.3901.002
----------------	-------------

## Product descriptions

Name	Tissue culture dish, (ØxH): 60 x 15 mm, Surface: Standard
Product number	83.3901.002
Product description	Tissue culture dish, (ØxH): 60 x 15 mm, with grid, material: PS, surface: Standard, for adherent cells, coding colour: red, TC Tested, 10 piece(s)/bag

## Material & colours

Closure material	Polystyrene (PS)
------------------	------------------

## Size

Diameter	60 mm
Height of product	15 mm

## Product characteristics

Growth area	22.1 cm <sup>2</sup>
Base shape	SUREGrip Design
Vessel type	Dish
Closure type	Slip-on lid
Cell type	adherent cells

## Purity & certification

Product category	no medical device   IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Piece(s) / inner box	10
Piece(s) / case	500
Piece(s) / pallet	12000
Height of case	380 mm
Width of case	332 mm
Depth of case	332 mm
Volume of case	0.0419 cbm
Weight of case	4.02 kg
EAN inner box	4038917338103
EAN case	4038917338097

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT