




Product specification

**Tissue culture dish, (ØxH): 100 x 20 mm, surface:
Standard**



Product information

| | |
|----------------|---------|
| Product number | 83.3902 |
|----------------|---------|

Product descriptions

| | |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------|
| Name | Tissue culture dish, (ØxH): 100 x 20 mm, Surface: Standard |
| Product number | 83.3902 |
| Product description | Tissue culture dish, (ØxH): 100 x 20 mm, material: PS, surface: Standard, for adherent cells, coding colour: red, TC Tested, 10 piece(s)/bag |

Material & colours

| | |
|------------------|------------------|
| Product material | Polystyrene (PS) |
| Closure material | Polystyrene (PS) |

Size

| | |
|-------------------|--------|
| Diameter | 100 mm |
| Height of product | 20 mm |

Product characteristics

| | |
|-------------|----------------------|
| Growth area | 58.8 cm ² |
| Base shape | SUREGrip Design |
| Vessel type | Dish |

Closure type

Slip-on lid

Cell type

adherent cells

Purity & certification

Product category

no medical device | IVD

Purity standard

TC Tested

Sterilisation

Gamma irradiation

Pyrogen-free

yes

Batched

yes

Packaging

Minimum order quantity

300

Piece(s) / inner box

10

Piece(s) / case

300

Piece(s) / pallet

4200

Height of case

460 mm

Width of case

313 mm

Depth of case

505 mm

Volume of case

0.0727 cbm

Weight of case

7.3701 kg

EAN inner box

4038917338196

EAN case

4038917338189

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT