




Product specification

Tissue culture dish, (ØxH): 100 x 20 mm, surface: Cell+



Product information

Product number	83.3902.300
----------------	-------------

Product descriptions

Name	Tissue culture dish, (ØxH): 100 x 20 mm, Surface: Cell+
Product number	83.3902.300
Product description	Tissue culture dish, (ØxH): 100 x 20 mm, material: PS, surface: Cell+, for challenging adherent cells, coding colour: yellow, TC Tested, 10 piece(s)/bag

Material & colours

Closure material	Polystyrene (PS)
------------------	------------------

Size

Diameter	100 mm
Height of product	20 mm

Product characteristics

Growth area	58.8 cm ²
Base shape	SUREGrip Design
Vessel type	Dish
Closure type	Slip-on lid
Cell type	challenging adherent cells

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	10
Piece(s) / case	300
Piece(s) / pallet	4200
Height of case	450 mm
Width of case	313 mm
Depth of case	505 mm
Volume of case	0.0711 cbm
Weight of case	7.48 kg
EAN inner box	4038917338257
EAN case	4038917338240

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT