



Product specification

Cell culture flask, T-25, surface: Cell+, 2-position screw cap



Product information

Product number	83.3910.300
----------------	-------------

Product descriptions

Name	Cell culture flask, T-25, Surface: Cell+, 2-position screw cap
Product number	83.3910.300
Product description	Cell culture flask, T-25, material: PS, surface: Cell+, for challenging adherent cells, coding colour: yellow, 2-position screw cap, TC Tested, 10 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Closure material	Polyethylene (PE)
Color of cap	yellow

Product characteristics

Growth area	25 cm ²
Vessel type	Flask
Cap	2-position screw cap
Cell type	challenging adherent cells

Purity & certification

Product category	no medical device IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	300
Piece(s) / inner box	10
Piece(s) / case	300
Piece(s) / pallet	6000
Height of case	340 mm
Width of case	408 mm
Depth of case	602 mm
Volume of case	0.0835 cbm
Weight of case	7.42 kg
EAN inner box	4038917337410
EAN case	4038917337403

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT