



## Product specification

### Cell culture flask, T-25, surface: Suspension, 2-position screw cap



## Product information

Product number	83.3910.500
----------------	-------------

## Product descriptions

Name	Cell culture flask, T-25, Surface: Suspension, 2-position screw cap
Product number	83.3910.500
Product description	Cell culture flask, T-25, material: PS, surface: Suspension, for suspension cells, coding colour: green, 2-position screw cap, TC Tested, 10 piece(s)/bag

## Material & colours

Closure material	Polyethylene (PE)
Color of cap	green

## Product characteristics

Growth area	24.7 cm <sup>2</sup>
Vessel type	Flask
Cap	2-position screw cap
Cell type	suspension cells

## Purity & certification

Product category	not regulated
------------------	---------------

CE certified	no CE mark
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Piece(s) / inner box	10
Piece(s) / case	300
Piece(s) / pallet	6000
Height of case	350 mm
Width of case	408 mm
Depth of case	595 mm
Volume of case	0.085 cbm
Weight of case	7.5 kg
EAN inner box	4038917337779
EAN case	4038917337762

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT