




Product specification

miniPERM®, Filling tube 5", for miniPERM® bioreactor



Product information

Product number	94.6077.138
----------------	-------------

Product descriptions

Name	miniPERM®, Filling tube 5", for miniPERM® bioreactor
Product number	94.6077.138
Product description	miniPERM®, Filling tube 5", for miniPERM® bioreactor, not autoclavable, sterile, 1 piece(s)/blister

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	5.9 mm
Width of product	5.9 mm
Length of product	5.9 mm
Height of product	127 mm

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile
Sterilisation	Gas sterilisation

Batched	yes
---------	-----

Packaging

Piece(s) / inner box	50
Piece(s) / case	300
Piece(s) / pallet	14400
Height of inner box	80 mm
Width of inner box	150 mm
Depth of inner box	200 mm
Height of case	260 mm
Width of case	268 mm
Depth of case	353 mm
Volume of case	0.0246 cbm
Weight of case	1.5501 kg
EAN inner box	4038917376488
EAN case	4038917255974

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT