



## Product specification

### Cell culture flask, T-75, surface: Standard, Filter cap



## Product information

Product number	83.3911.002
----------------	-------------

## Product descriptions

Name	Cell culture flask, T-75, Surface: Standard, Filter cap
Product number	83.3911.002
Product description	Cell culture flask, T-75, material: PS, surface: Standard, for adherent cells, coding colour: red, Filter cap, TC Tested, 5 piece(s)/bag

## Material & colours

Product material	Polystyrene (PS)
Closure material	Polyethylene (PE)
Color of cap	red

## Product characteristics

Growth area	75 cm <sup>2</sup>
Vessel type	Flask
Cap	Filter cap
Cell type	adherent cells

## Purity & certification

Product category	no medical device   IVD
------------------	-------------------------

Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

# Packaging

Minimum order quantity	100
Piece(s) / case	100
Piece(s) / pallet	1600
Height of case	435 mm
Width of case	402 mm
Depth of case	445 mm
Volume of case	0.0778 cbm
Weight of case	6.78 kg
EAN inner box	4038917337687
EAN case	4038917366052

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**